# BILIRRUBINA TOTAL

## MÉTODO DE JENDRASSIK – GROF

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Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

# TOTAL BILIRUBIN

## JENDRASSIK – GROF METHOD

For "in vitro" determination of total bilirubin in serum or plasma

### PRINCIPLE

Total bilirubin is determined by reaction with diazotized sulphanilic acid, in the presence of caffeine, with the final production of an azopigment. The same reaction, but in the absence of caffeine, is used to measure direct bilirubin.

### DIAGNOSTIC USE

Bilirubin is a metabolic hemoglobin product which is transported to the liver, where there is conjugated with glucuronic acid and it is secreted through bile. This conjugated bilirubin is called direct bilirubin, while unconjugated bilirubin is called indirect bilirubin. The total serum bilirubin is direct bilirubin plus indirect bilirubin.

Direct bilirubin increases in biliary obstruction, cirrhosis and hepatitis. High level of indirect or total bilirubin may indicate hemolytic diseases, physiological neonatal jaundice, Gilbert's disease or fructose intolerance.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

**Kit 1 x 250 mL. (Ref. 99 27 14).** Contents:

- A. 1 x 50 mL Sulphanilic acid. Ref. 99 90 11
- B. 2 x 100 mL Caffeine. Ref. 99 23 91
- C. 1 x 2 mL Sodium nitrite. Ref. 99 90 94

### WORKING REAGENT PREPARATION

Preparation of the diazo reagent: Mix 1 part of solution **C** with 50 parts of solution **A**.

### REAGENTS COMPOSITION

The concentrations in the reagent solutions are:

- A.** Sulphanilic acid 6 mM  
HCl 0.17 M
- B.** Caffeine 0.25 M  
Sodium benzoate 0.50 M
- C.** Sodium nitrite 0.4 M

### STORAGE AND STABILITY

The components of the kit, when stored at room temperature ( $\leq 25^{\circ}\text{C}$ ), will remain stable until the expiration date stated on the label.

Diazo reagent is stable for 15 days at room temperature ( $\leq 25^{\circ}\text{C}$ ).

### Signs of reagent deterioration:

Presence of particles or turbidity in the reagent. Working reagent blank  $>0.100$ .

### ADDITIONAL EQUIPMENT

General laboratory equipment.

Spectrophotometer or photometer. Cuvette: 1 cm light-path

### SAMPLE

Serum or plasma free from haemolysis. Serum bilirubin will decrease 50% in one hour if kept at room temperature ( $\leq 25^{\circ}\text{C}$ ), and at direct sunlight. Samples kept in the dark at  $2-8^{\circ}\text{C}$  will remain stable for up to 3 days.

### CAUTION

The reagents contains Sodium azide at 0.09%. Handle with care.

The safety statements are on the label. We advise to read MSDS before reagent handling.

Waste products must be handled as per local regulations.

### AUTOANALYZERS

Adaptations to different autoanalyzers are available on request.

### PROCEDURE

Technique	BL S (mL)	SA (mL)
Reagent B (Caffeine)	1.0	1.0
Diazo reagent	--	0.1
Saline	0.1	--
Sample	0.1	0.1

Mix and let stand at room tem ( $20 - 25^{\circ}\text{C}$ ). **Read exactly at 10 minutes.**

### Reading

Wavelength: 546 nm

Blank:  $\text{H}_2\text{O}$

Cuvette: 1 cm light-path

### CALCULATIONS

$\Delta \text{Abs} = \text{Abs SA} - \text{Abs BL S}$

Where:

Abs SA: Sample Absorption

Abs BLS: Sample blank Absorption

$\Delta \text{Abs} \times 12.8 = \text{mg total bilirubin} / \text{dL}$

### SI Units

$\text{mg} / \text{dL} \times 17.1 = \mu\text{mol} / \text{L}$

### REFERENCE VALUES

Adults: up to 1 mg/dL. total bilirubin

Neonate:	Premature	Full term
< 1 day	1.0 - 8.0 mg/dL	2.0 - 6.0 mg/dL
< 2 days	6.0 - 12.0 mg/dL	6.0 - 10.0 mg/dL
3-5 days	10.0 - 14.0 mg/dL	4.0 - 8.0 mg/dL

The stated values are for guidance. Each particular laboratory should establish its own normal range.

### PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using an automated analyzer.

Sensitivity, as detection limit: 0.1 mg/dL

Linearity: Up to 25 mg/dL. Concentrations higher than 25 mg/dL may produce an absorbance value ( $\Delta \text{Abs} > 1.500$ ) that can hardly be read in non-digital photometers. It is recommended in such a case, to carry out a 1/10 dilution of the sample with saline, (NaCl 0.9%) and then to multiply the final result by 10.

Accuracy: 97.8%

Repetitivity, as CV%: 1.3%

Reproducibility, as CV%: 1.7%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the comparison experiments are available on request.

### INTERFERENCES

Avoid hemolysed samples.

### QUALITY CONTROL

Control serum, Seriscann Normal (Ref. 99 41 48) and Seriscann Abnormal (Ref. 99 46 85) should be included in each test series. Each particular laboratory should establish its own control program.

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