# APOLIPOPROTEINA B

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Richitit, RF. et. al., *Clin. Chem.*, 1991, 37, 1676.

Young, D.S.; Effects of Preanalytical Variables on Clinical Laboratory Tests. 3rd Ed. AACCPress 2000.

Young, D.S.; Effects of drugs on Clinical Laboratory Tests. 4th Ed. AACCPress. 2000.

# APOLIPOPROTEIN B

## QUANTITATIVE DETERMINATION BY IMMUNOTURBIDIMETRY

For "in vitro" diagnostic

### PRINCIPLE

The reagent allows quantifying the level of human apolipoprotein B present in the sample by comparing the turbidimetric response produced with that obtained from the standard curve of known concentrations of human apolipoprotein B.

The reagent is an antiserum anti-human apolipoprotein B which reacts with the apolipoprotein B of the serum giving protein aggregates. This aggregation process produces an increase in the Abs. of the system.

### DIAGNOSTIC USE

The apolipoprotein B (APO B) is the most abundant structural protein in LDL, VLDL lipoprotein and chylomicrons. The main function is to transport lipoproteins rich in triglycerides from the liver to body tissues.

Since changes on APO A1 and APO B concentrations are similar to the ones in HDL and LDL on patients with coronary heart diseases, the measurement of apolipoproteins can be applied to the estimation of coronary risk.

Increase in APO B concentration could be because a disorder on lipoprotein metabolism as hyperapobetalipoproteinemia, in which the concentration of APO B increase and LDL-cholesterol levels are normal.

APO B concentration decreases in abetalipoproteinemia, in which the secretion of intestinal and hepatic lipoproteins rich in triglycerides is blocked.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

**Kit (Ref. 99 34 06).** Contents:

- A. 1 x 16.0 mL Buffer solution. Ref. 99 04 59  
B. 1 x 4.0 mL anti-human APO B serum Ref. 99 04 83

### Apolipoprotein Calibrator.

Ref. 99 36 90

1 x 1 mL Pool of human plasma samples calibrated in apolipoprotein A1 and apolipoprotein B.

The concentration of each component is stated on the vial label.

### Apolipoprotein Control

Ref. 99 36 92

1 x 1 mL Samples positive for apolipoprotein A1 and apolipoprotein B.

The concentration of each component is stated on the vial label.

### REAGENT PREPARATION

Reagents are ready to use.

### REAGENT COMPOSITION

- A. Sodium phosphate/saline buffer pH 7.4 with stabilizers and preservatives.  
B. Goat antiserum anti-human APO B in buffered saline pH 7.4 with stabilizers and preservatives.

### STORAGE AND STABILITY

The antiserum reagent, the buffer solution and the Calibrator and Control are stable up to the expiry date indicated on the label when kept at 2-8°C.

Do not freeze antisera or buffer solution.

### Signs of reagent deterioration:

Presence of particles in the reagent. Quality control values outside acceptance range.

### ADDITIONAL EQUIPMENT

General laboratory equipment.  
Spectrophotometer or photometer thermostatable at 37°C with a 340nm filter.

### CAUTION

Human plasmas used in the control preparation have been found negative in the reaction to HBsAg and HIV I/II. However they should always be handled with care. All reagents contain sodium azide 0.09% as preservative. It is recommended to read SDS before reagent manipulation.

Waste products must be handled as per local regulations.

### SAMPLE

Recent serum or plasma (heparin or EDTA plasma) kept at 2-8°C for no longer than 48h. If the test is to be performed later, it is recommended to freeze the sample. Avoid successive freezing and thawing. Discard hemolyzed or contaminated samples.

### QUALITY CONTROL

Calibrator or positive control should be included in each run. Each particular laboratory should establish its own control program.

### PROCEDURE

#### Technique

To be used in clinical chemistry analyzers.

1. Shake gently the reagents before use. The sample shall be processed without dilution.

Conditions:

Volume of buffer solution (A):	240 $\mu$ L
Volume of antibody solution (B):	60 $\mu$ L
Volume of sample:	3 $\mu$ L
Reaction temperature:	37°C

2. Mix well; start the stopwatch and read the Abs (Absi). Repeat the reading after 2min of reaction (Absf) at 340nm

### CALCULATIONS

Determine the  $\Delta$ Abs (Absorbance increase) as (Absf – Absi) for each sample, standard or control. Apolipoprotein B concentrations are obtained by the interpolation of the  $\Delta$ Abs value in the calibration curve.

#### Calibration curve

Prepare the curve by serially diluting the calibrator with saline (apolipoprotein calibrator Ref. 99 36 90). To obtain the concentration of each dilution, multiply the calibrator concentration by the corresponding factor indicated in the table:

	0	1	2	3	4	5
CAL ( $\mu$ L)	--	25	50	100	200	400
Saline sol. ( $\mu$ L)	400	375	350	300	200	--
Factor	0	0.0625	0.125	0.25	0.5	1.0

The method described above is only indicative, as it should be specifically adapted to each autoanalyzer.

### REFERENCES VALUES

59-168 mg/dL

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

### PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using an autoanalyzer.

Sensitivity, as detection limit: 7.8 mg/dL

Reaction range: 300 mg/dL. Samples that give higher concentration should be diluted in saline NaCl 0.9% (1+1) and the final result has to be multiplied by 2.

Accuracy: 95.1%

Repetitivity, as CV%: 1.69%

Reproducibility, as CV%: 7.19%

No prozone effect was observed up to 3000 mg/dL

Trueness: Results obtained with this reagent did not show systematic differences when compared with the reference reagent.

Details of the comparison experiments are available on request.

### INTERFERENCES

No interference was observed by hemoglobin (125 mg/dL), bilirubin (5 mg/dL) and RF (600 UI/mL). Lipemic sera and other substances could interfere (See references).

Samples that give an Abs. higher than that obtained for the last point in the standard curve should be diluted in saline (NaCl 0.9%) and the determination repeated.

### AUTOANALYZERS

QCA has adaptation to different autoanalyzers on request.

### REFERENCES

Tiezt, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012).

Richitie, RF. et. al., *Clin. Chem.*, 1991, 37, 1676.

Young, D.S.; Effects of Preanalytical Variables on Clinical Laboratory Tests. 3rd Ed. AACCPress 2000.

Young, D.S.; Effects of drugs on Clinical Laboratory Tests. 4th Ed. AACCPress. 2000.

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