


A.F. Genital System

System for detection, count and susceptibility testing of urogenital mycoplasmas and pathogenic microorganisms

DIRECT INOCULATION FROM THE SPECIMEN SUSPENSION. **NO ENRICHMENT BROTH NEEDED**



RESULTS IN 24-48 HOURS INCUBATION AT 36±1 °C INEXPENSIVE, EASY TO USE, INSTRUMENTS FREE

Urogenital Mycoplasmas

COUNT

<104 CFU/mL



Low growth of mycoplasmas is shown by colour change of well: from yellow to red



Moderate growth of mycoplasmas is shown by colour change of wells: I-GR+, 2-GR++ from yellow to red



High growth of mycoplasmas is shown by colour change of wells: I-GR+, 2-GR++, 3-GR+++ from yellow to red

DETECTION

Mycoplasma hominis



The presence of Mycoplasma hominis is shown by colour change of well: 4-ADC from yellow to red

Ureaplasma urealyticum



The presence of Ureaplasma urealyticum is shown by colour change of well: 5-UR from yellow to red

Mycoplasma hominis and Ureaplasma urealyticum



The presence of Micoplasma hominis and Ureaplasma urealyticum is shown by colour change of wells: 4-ADC and 5-UR from yellow to red

SUSCEPTIBILITY TESTING





















The **Sensitivity** to antibiotics is shown by the colour of wells: vellow



The Intermediate sensitivity to antibiotic is shown by colour change of wells: from yellow to orange



The **Resistance** to antibiotic is shown by colour change of wells: from yellow to red

Presumptive Identification

Trichomonas vaginalis / Yeasts



Observe the presence of mobile ciliated trophozoites, chlamydospores and hyphae 6-TRYE under microscopic observation (40x).





Trichomonas vaginalis trophozoite.



Candida albicans chlamydospores and

Escherichia coli



The presence of **Escherichia coli** is shown by colour change of well: 16 - ESC from red to blue



The presence of **Proteus spp.** is shown by colour change of well: 17 - PRO from yellow to brown 17 - PRO



The presence of Pseudomonas spp.

is shown by colour change of well: 18 - PSE from yellow to green



Gardnerella vaginalis



The presence of Gardnerella vaginalis is shown by colour change of well: 19 - GAR from red to yellow



Staphylococcus aureus The presence of Staphylococcus aureus is shown by colour change of well:

20 - STF from yellow to black



Enterococcus faecalis The presence of **Enterococcus faecalis** is shown by colour change of well:



Neisseria gonorrhoeae



The presence of **Neisseria** gonorrhoeae is shown by colour change of Oxidase test stick: from white to blue



Streptococcus agalactiae

The presence of Streptococcus agalactiae is shown by colour change of well:

23 - STG from yellow to green





The presence of Candida spp. is shown by colour change of well: 24 - CAN from green to yellow 24 - CAN

21 - STR from yellow to black





A.F. Genital System


TEST PROCEDURE

Code 74156

20 test


EXAMPLE OF RESULTS





Sample of vaginal swab with presence of **Gardnerella vaginalis** (19-GAR) **Enterococcus faecalis** (21-STR) **Streptococcus agalactiae** (23-STG) **Candida spp.** (24-CAN)

(2)



Sample of vaginal swab with presence of **Mycoplasma hominis** (4-ADC) **Ureaplasma urealyticum** (5-UR) with growth +++ (>10⁵ CFU/mL) **E. coli** (16-ESC) and **Gardnerella vaginalis** (19-GAR)

Susceptibility testing of mycoplasmas: II-E, I2-CLA, I4-IOS, I5-CD = Resistant



Immerse the swab (after obtaining the clinical material) in the vial of physiological solution and wait 5 minutes. Carefully squeeze the swab against the vial wall.



Cover wells I to 5, 7 to 15, 19 and 24 with I drop of Vaseline Oil. Cover the system with the lid provided and incubate in thermostat at 36 ± 1 °C for 18-24 hours.



Take a drop of liquid from well **22-NES**, deposit in on an Oxidase Test Stick and observe the immediate development (10 seconds roughly) of a blue colour (positive oxidase test)



Transfer 0.2 mL of clinical sample suspension into each well of the system.



Take a drop of liquid from well **6-TR/YE**, deposit it on a glass slide, place a cover slip on top and examine under the microscope (40x) for the presence of *Trichomonas vaginalis* and *Candida* spp.



Read and interpret by automatic reader or manually with the test results form.



 ϵ

