


URIN SYSTEM Plus

System for count, presumptive identification and susceptibility testing of microorganisms from urinary infections

DIRECT INOCULATION FROM URINE SAMPLE



RESULTS IN 24-48 HOURS INCUBATION AT 36±1 °C **INEXPENSIVE, EASY TO USE, INSTRUMENTS FREE**

Total Bacterial Count

NO microbial growth



I-GR+, 2-GR+: violet No colour change

Significant bacteriuria



I-GR+: from violet to yellow Significant bacteriuria (105 < CFU/mL < 106)

Marked bacteriuria



I-GR+, 2-GR++: from violet to yellow Marked bacteriuria (CFU/mL>106)

Presumptive Identification

Escherichia coli



The presence of **Escherichia coli** is shown by colour change of well: 3 - ESC from red to blue




Proteus spp. The presence of Proteus spp.

is shown by colour change of well:

4 - PRO from yellow to brown



Pseudomonas spp.





KES group



The presence of is shown by colour change of well: 6 - KES from violet to yellow



Enterococcus spp.

The presence of Enterococcus spp. is shown by colour change of well: 7 - STR 7 - STR from yellow to black



Staphylococcus aureus







The presence of Candida spp. is shown by colour change of well:



Observe the precence of chlamydospores and hyphae under microscopic observation (40x)



Candida albicans chlamydospores and hyphae

Susceptibility Testing





The **Sensitivity** to antibiotics is shown by colour of wells: blue



The Intermediate sensitivity to antibiotics is shown by colour change of wells: from blue to grey



The **Resistance** to antibiotics is shown by colour change of wells: from blue to yellow



The Microbial Growth Control is shown by colour change of well: 24-C from blue to yellow



URIN SYSTEM Plus


TEST PROCEDURE

Code 74160

20 test

EXAMPLE OF RESULTS






Sample of urine with Marked bacteriuria (I-GR+, 2-GR++) and presence of Proteus spp. (4-PRO) and Enterococcus spp. (7-STR)

Susceptibility testing: 10-AK, 11-CN, 12-TOB, 16-CTX, 17-CAZ, 18-AMS, 19-NA, 23-SXT = **Resistance**

(2)



Sample of urine with **Marked bacteriuria** (I-GR+, 2-GR++) and presence of E. coli (3-ESC)

Susceptibility testing: 23-SXT = **Resistance**




add 0.5 mL of urine to the vial of Physiological Solution.



- Transfer 0.2 mL of **Suspension A** to each well from 1 to 9. (Total bacterial count and identification)
- Transfer 0.2 mL of **Suspension B** to each well from 10 to 24 (Susceptibility testing)



to the vial of Suspension Medium.



- Incubate in thermostat at 36±1 °C for 18-24 hours.
 At the end of incubation, watch for colour change of the wells and interpret the results.



Take a drop of liquid from well **9-CAN**, deposit it on a glass slide, place a cover slip on top and examine under the microscope (40x) for the presence of

Candida spp. (chlamydospores and hyphae).



Read and interpret by automatic reader or manually with the test results form.



