

CONTACT SLIDE



FLEXIBLE CONTACT SLIDE

FOR BACTERIOLOGICAL
SCREENING OF SURFACES

ISO 18593

DESCRIPTION	PACKAGING	REF.
CONTACT SLIDE 1 (PCA+TTC + NEUTRALIZING / V.R.B.G. AGAR + NEUTRALIZING) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526
CONTACT SLIDE 2 (PCA+TTC + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527
CONTACT SLIDE 3 (PCA+TTC + NEUTRALIZING / VOGEL JOHNSON AGAR) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530
CONTACT SLIDE 4 (CETRIMIDE AGAR + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) <i>Pseudomonas</i> , yeasts and moulds.	20 slides 120 slides	525282 53528
CONTACT SLIDE 5 (V.R.B.G. AGAR + NEUTRALIZING / BILE ESCULIN AGAR + NEUTRALIZING) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531
CONTACT SLIDE 8 (PCA+TTC + NEUTRALIZING / V.R.B.L. AGAR + NEUTRALIZING) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562
CONTACT SLIDE 9 (PCA+TTC + NEUTRALIZING / PCA+TTC + NEUTRALIZING) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532



microbiology products

FS0616
Rev. 3 - 13/11/2010